

# MULTI-THERM™

COOL - HEAT - SHAKE



- Two models for temperature control from 0° to 100°C
- Exchangeable blocks, for tube sizes 0.2 to 50ml
- Precision contoured wells, for uniform thermal transfer

Temp. Range (HC): Ambient -20 to 100°C  
 Temp. Range (H): Ambient +5 to 100°C  
 Speed Range: 200 to 1500 rpm  
 Shaking Orbit: 3mm  
 Temp. Accuracy: +/- 1.0°C  
 Dimensions (LxWxH): 8.25 x 11.5 x 8.3 in.

Wide Assortment of Block Options



**Ordering Information:**

|                                 |                                      |               |
|---------------------------------|--------------------------------------|---------------|
| H5000-HC                        | <b>MULTI-THERM</b> Shaker            | <b>\$2495</b> |
| H5000-H                         | <b>MULTI-THERM</b> Shaker, heat only | <b>\$1695</b> |
| Please Order Blocks Separately: |                                      |               |
| H5000-02                        | 96 x 0.2ml or one PCR plate          | \$295.00      |
| H5000-05                        | 54 x 0.5ml                           | \$249.00      |
| H5000-CMB                       | Combi, 15x0.5ml & 20x1.5ml           | \$249.00      |
| H5000-15                        | 35 x 1.5ml                           | \$249.00      |
| H5000-20                        | 35 x 2.0ml                           | \$249.00      |
| H5000-12                        | 24 x 12mm*                           | \$349.00      |
| H5000-150                       | 12 x 15ml*                           | \$349.00      |
| H5000-500                       | 6 x 50ml*                            | \$349.00      |
| H5000-MP                        | 1 x Microplate*                      | \$349.00      |

\* Limited maximum speed of 750rpm, minimum temperature of ambient -10C

Order A Multi-Therm or Incu-Shaker Mini before 11/30/11 and receive one of the items below . . . FREE!

**Special Offer:**

Promotional items must be included on order with H5000-H, H5000-HC or H1000-M



**FREE!**

Cat No. R1000  
Benchtop Cooler



**FREE!**

Cat No. BV1000  
Vortex Mixer



**FREE!**

Cat No. C1008-B  
Mini Centrifuge



**Personal Sized**

Temp. Range: Ambient +5 to 60°C  
 Shaking Speed: 30 to 300rpm  
 Shaking Orbit: 19mm (3/4 in.)  
 Dimensions (LxWxH): 15.7 x 11 x 10.5 in.

# INCUSHAKER™ Mini

- Exceptional temperature uniformity
- Extremely compact, just 11" wide
- Exchange flask clamps in one second with new MAGIC CLAMP™ system



MAGIC CLAMP™ Instant Clamp Attachment

**Ordering Information:**

|               |  |           |
|---------------|--|-----------|
| H1000-M       | <b>INCUSHAKER™</b> Mini Shaking Incubator with non-slip rubber mat                     | \$2495.00 |
| H1000-MR      | <b>MAGIC CLAMP</b> universal platform for flasks & tube racks (clamps sold separately) | \$209.00  |
| H1000-MR-500  | Clamp for Magic Clamp platform, 500ml (max. 5)   | \$27.00   |
| H1000-MR-250  | Clamp for Magic Clamp platform, 250ml (max. 8)   | \$26.00   |
| H1000-MR-125  | Clamp for Magic Clamp platform, 125ml (max. 12)  | \$24.00   |
| H1000-MR-1550 | Tube rack for Magic Clamp platform, 12x50 and 21x15 ml conical tubes (max. 1)          | \$129.00  |

Additional accessories are also available. Visit [www.benchmarkscientific.com](http://www.benchmarkscientific.com) for a complete listing.