

hybex™

Erlenmeyer Flasks



hybex™ Erlenmeyer Flasks are manufactured from high quality, **borosilicate 3.3 glass**. A low coefficient of thermal expansion (3.3), provides temperature resistance that is far superior to that of standard lab glass (or soda lime glass). Hybex glass is designed to withstand high temperatures, making them ideal for shaking applications, media preparation and long-term storage.

The standard Erlenmeyer flasks feature a reinforced rim along with permanent, white enamel graduations for volume measurement and a white identification box for marking or coding.

For storage of media and cultures, the capped flasks include an autoclavable polypropylene cap & drip-ring. Optional cap colors are available, as well as PBT material caps for temperature applications up to 180°C. (Threads are standard GL-45 for compatibility with a wide range of cap types from various manufacturers, including phenolic caps, vented caps, etc.)

- Available w/GL45 screw cap (and drip rings)
- Borosilicate 3.3 hybex glass
- Autoclavable, reuseable & breakage resistant
- Permanent white enamel marking area
- Ideal for shaking, media prep or storage



Available in standard (non-threaded)

Technical Data

Material:	hybex™ borosilicate glass (3.3)
Cap Thread Size:	GL45*
Sterilization:	Autoclave or dry heat
Cap Material:	Polypropylene Autoclavable (up to 140°C)
Included Cap:	Blue

Ordering Information

Erlenmeyer Flasks, Borosilicate 3.3 Glass (standard, no caps or threads)

B3110-50	50ml Erlenmeyer flask, pack of 12
B3110-125	125ml Erlenmeyer flask, pack of 12
B3110-250	250ml Erlenmeyer flask, pack of 12
B3110-500	500ml Erlenmeyer flask, pack of 6
B3110-1000	1000ml Erlenmeyer flask, pack of 6
B3110-2000	2000ml Erlenmeyer flask, 1 flask

Erlenmeyer Flasks, Borosilicate 3.3 Glass with GL45 Caps

B3220-250	250ml Erlenmeyer flask with GL45 cap, pack of 12
B3220-500	500ml Erlenmeyer flask with GL45 cap, pack of 6
B3220-1000	1000ml Erlenmeyer flask with GL45 cap, pack of 6
B3220-2000	2000ml Erlenmeyer flask with GL45 cap, 1 flask
B3000-HTC	High temperature red caps, up to 180°C, PBT (GL45), 10/pk
B3000-CAP*	Replacement cap, standard blue (GL45), 10/pk.

*Please add (-G) for green caps, (-P) for purple caps, (-Y) for yellow caps, (-W) for white caps or (-AS) for assorted caps (2 of ea. color)

Benchmark
Scientific



Benchmark Scientific Inc
PO Box 709 Edison, NJ 08818
Phone: 908-769-5555 Fax: 908-222-1864
Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com