

# Sprint™ 6H Portable

by Benchmark



6 position swing out rotor



## Clinical testing on the go!

The Sprint™ 6H Portable kit allows for mobile processing of clinical applications. The package includes a 6 position swing out rotor for blood tubes, 12V voltage adapter, and adjustable leveling stands to match various car seat angles and ensure balanced centrifugation.

Speed:	3,400rpm / 1,600 xg
Capacity:	6 x 10ml 6 X 5-7ml (adapters required)
Timer:	1 to 99 mins.
Dimensions:	12.8 x 11.5 x 9 in. (W x D x H) (32.5 x 29.3 x 23 cm)
Weight:	17.6 lbs / 8 kg
Electrical:	120 or 230V, 50-60 Hz,
Warranty:	2 years



## Ordering information

C5000-6H-PT Sprint™ 6H portable kit  
Includes: 12V power supply, adjustable stands

**Benchmark** Scientific

PO Box 709 Edison, NJ 08817  
TEL: 908-769-5555 Fax: 732-313-7007  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)