

Now is the time to stock up on our most popular QuickSilver Buffers  
With our **BUY ONE GET ONE FREE\*** offer



Our exclusive quick-dissolving formulation and premeasured pouches offer the utmost in convenience.

Mix one pouch with 1 liter of distilled water and you have fresh 1X buffer that's ready to use.


- Perfectly formulated, fresh buffer ready when you need it
- No more measuring, weighing and mixing components
- Guaranteed quality and consistency



Convenient, premeasured packs



Just add pouch contents to water and stir to dissolve



Supplied with preprinted bottle labels



**Buy 1**

**Get 1 FREE!**

**Buffer  
PROMO**

**\*Buy one case at regular price and choose an additional case of equal or lesser value for FREE.**

Use promo code BB2021. Offer valid through 10/31/21

Cat#	Description	List Price
EB1000	QuickSilver TAE Buffer Powder, 100 pouches	\$204.00
EB1001	QuickSilver TBE Buffer Powder, 50 pouches	\$165.00
EB1002	QuickSilver Fast Running Buffer Powder, 100 pouches	\$187.00
EB1100	QuickSilver Tris-Glycine-SDS Buffer Powder, 50 pouches	\$145.00
EB1101	QuickSilver Tris-MES-SDS Buffer Powder, 50 pouches	\$295.00
EB1102	QuickSilver Tris-MOPS-SDS Buffer Powder, 50 pouches	\$295.00
EB1202	QuickSilver TBS Buffer Powder, 100 pouches	\$180.00
EB1203	QuickSilver TBS-T Buffer Powder, 100 pouches	\$180.00

\*Each pouch makes 1 liter of 1X buffer




**FREE BUFFER** (up to \$204 value) with the purchase of any of our Mini Electrophoresis Systems:



## myGel™ Mini Electrophoresis System

- Complete system, includes gel box, power supply and casting sets
- Quick and easy gel casting, no tape, clamps or spacers required
- Compact design, footprint: 7 x 5.1 in.

Cat#	Description	List Price
E1101	myGel™ Mini Electrophoresis System	\$545.00



## myGel InstaView™ Electrophoresis System

includes direct-connect mini power supply

- View DNA migration and separation as it happens
- Blue light illumination – won't damage DNA
- Capture high resolution images without removing the gel from the tank

Cat#	Description	List Price
E2101	myGel InstaView™ Complete Electrophoresis System	\$875.00




Smartphone imaging enclosure included

## myGel InstaView™ Electrophoresis Tank System

connects to any standard power supply

- View DNA migration and separation as it happens
- Blue light illumination – won't damage DNA
- Capture high resolution images without removing the gel from the tank



E1200 shown connected to myVolt Touch Power Supply, sold separately

Cat#	Description	List Price
E2100	myGel InstaView™ Gel Tank with Blue LED Illuminator	\$670.00
E2301	myVolt™ Touch Power Supply	\$540.00



Smartphone imaging enclosure included



Choose one case of our Quick Silver TAE, TBE or Fast Running Buffer with the purchase of E1101, E1201, or E1200 system above.

Cat#	Description	List Price
EB1000	QuickSilver TAE Buffer Powder, 100 pouches	\$204.00
EB1001	QuickSilver TBE Buffer Powder, 50 pouches	\$165.00
EB1002	QuickSilver Fast Running Buffer Powder, 100 pouches	\$187.00

\*Each pouch makes 1 liter of 1X buffer