## BenchMixer Benchmark Vortex Mixer



- **Powerful Motor for INSTANT Vortexing**
- Q-Drive<sup>™</sup> Dynamic Balance System
- **Touch or Continuous Operation**
- Variable Speed Control from 200 to 3200rpm











## Technical Data:

Speed Range: Operating Modes: Operating Temp. Range: **Dimensions:** 

Weight: Electrical:

Warranty:

200 to 3200 rpm **Touch or continuous** +4 C to +65 C 13 x 16 x 17 cm 5 x 6.3 x 6.7 in. 3.8 kg / 8.4 lb 115V AC, 60 Hz, 1.5A

230V AC, 50 Hz, 0.75A

2 years

The BenchMixer<sup>™</sup> sets the new standard in vortex mixing. With advanced Q-Drive <sup>™</sup> technology, the BenchMixer<sup>™</sup> provides smooth instant vortexing of even the largest sample sizes, including full 50ml tubes. The unique counter balance system creates maximum vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation.

Other features include: "no-walk" suction feet, cold room/incubator compatibility and spill-proof electronic controls.

## Ordering Information

BenchMixer Vortex Mixer BV1000\*

## **Optional Accessories:**

BV1000-FLAT Optional 3 inch flat head

BV1000-COMBO Optional head for tubes and microplate.

66 tubes (38 x 1.5ml and 28 x 0.5ml) or 1 microplate

BV1000-H15 Optional head for 12 x 1.5/2.0ml, horizontal

BV1000-H150 Optional head for 4 x 15ml, horizontal BV1000-H500 Optional head for 2 x 50ml, horizontal

\* To order a product in 230V, please add -230 to item number.



**Benchmark Scientific Inc.** Edison, NJ 08817 Phone: 908-769-5555 Fax: 908-222-1864

Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com