

Z446-K Refrigerated Universal Centrifuge

- Max. swing-out capacity:
28 x 50ml, 68 x 15ml, and 20 x microplates
- Temp Range -20°C to +40°C
- Bundle options for **tubes** or **plates!**
- Truly universal, 14 rotor options

As low as \$13,500



BUNDLE ORDERING INFORMATION*

		Reg.	Bundle
Z446-K-BND	Z446-K w/ Swing out rotor for 15/50ml	\$15,001	\$14,000
Z446-K-MPBND	Z446-K w/ Swing out rotor for microplates	\$15,045	\$13,500

Z216-MK Refrigerated High Speed Microcentrifuge

As low as \$5,900



- Temp Range -20°C to +40°C
- High speed, up to 21,631xg
- 2 innovative rotor bundle options!



BUNDLE ORDERING INFORMATION*

		Reg.	Bundle
Z216-K-BND44	Z216-MK Refrigerated w/ 44 place rotor	\$7,012	\$5,900
Z216-K-CMB	Z216-MK Refrigerated w/ COMBI-Rotor™	\$7,714	\$6,400

Ambient versions also available



Also available:
MC-24™ Touch
High Speed Microcentrifuge

- Combi-rotor for tubes and strips
24 x 1.5/2.0ml + 2 x PCR Strips
- Up to 16,800 xg!

\$2,375

BucketBuddy™

Magnetic Bucket and Tube Rack

Magnetic Modular

Accepts: Hermle inserts/buckets up to 750ml (150 mm x 145 mm area)

13 mm dia.
17 mm dia.
30 mm dia.

Part# Z326-BKR

*Additional bundle models available upon request. for 230V voltage, add "-E"

Offers valid until Dec. 31st 2022

FREE With ANY Hermle Universal Bundles
(Z306, Z326, Z366, Z446, Z32-HK and Z36-HK Series)



Unbeatable value!
\$4,500

Compact, only 14in. wide

Z287-A Micro Bundle

- Bundle includes 2 rotors
18 x 1.5 / 2.0 ml and 2 x microplate
- Up to 14,000 rpm and 16,058 xg
- Indicator Light, easily view cycle status
- 11 rotor options!



BUNDLE ORDERING INFORMATION*

Z287-BND-MB Z287-A with 2-place MP rotor and 18 x 1.5/2.0ml rotor

Reg. ~~\$5,254~~ **Bundle \$4,500**



As low as \$14,000

Z36-HK SUPER Speed Centrifuge

- Speeds up to 65,395 xg!
- Max. swing-out capacity: 12 x 50ml or 28 x 15ml
- 3 bundle options to meet any need!
- Refrigerated, -20°C to 40°C
- Includes FREE **BenchPACK!**™ lab starter pack! (**\$1,500 Value!**)



BUNDLE ORDERING INFORMATION*

	Reg.	Bundle
Z36-BND250 Z36-HK Centrifuge with 6 x 250ml rotor	\$19,718	\$15,800
Z36-BND50 Z36-HK Centrifuge with 6 x 50ml rotor	\$17,863	\$14,600
Z36-BND20 Z36-HK Centrifuge with 12 x 1.5/2.0ml rotor	\$16,967	\$14,000

Offers valid until Dec. 31st 2022

Additional Bundles Available



Z496-K Floor Standing*
w/ swing out rotor

\$26,000

*230V Only

SIEVA III Filtration Centrifuge
w/ perforated or non perforated basket

As low as \$22,000

Z206-A Compact Centrifuge
w/ 6 x 50ml rotor

\$2,600

Z306 Compact Universal
w/ swing out or microplate options

\$5,400

*Additional bundles available upon request.
To order 230V voltage, add "-E" to part number.