## StripSpin I2D Digital Mini Centrifuge



The StripSpin 12D from Benchmark is the only mini centrifuge designed for the spin down of 12-position PCR strips. The centrifuge features digital speed control up to 5,500 rpm, time control up to 99 min., and a safety interlock to stop the centrifuge if the lid is opened during operation.

With built in rcf conversion, speed can be set/displayed in rpm or rcf. The StripSpin 12D is the perfect solution for centrifugation of both 8 and 12-position PCR strips (or 48 x 0.2ml tubes) and includes a 2 year warranty.

Also available...

Includes COMBI-Roto for tubes and strips

Scient

FUGE 1294

Benchma

Item: C2212

Unique rotor accepts <u>12-position</u> PCR strips

- + Digital control up to 5,500 rpm
- "Mini" sized, less than 8 inches wide
- **H** Timer up to 99 mins.



## Holds up to 4 x PCR strips or 48 x 0.2ml tubes

Technical	<u>Data</u>	
Speed: Max. g-force: Capacity:	5,500rpm 1,690 xg 4 x 12 position PCR strips	
Dimensions:	4 x 8 position PCR strips 48 x 0.2ml PCR tubes 7.6 x 9 x 4.7 in	
(WxDxH) Weight:	19.4 x 23 x 12 cm 4.4 lbs / 2 kg	annun
Electrical: Warranty:	100 to 240V, 50-60Hz 2 Years	um.
	a fille	

## Ordering Information

C2248\*

*StripSpin* Mini Centrifuge, with rotor for up to 4 x PCR strips (8 or 12 position) or 48 x 0.2ml tubes

\*Includes US plug. For EU plug, add (-E) to the item number.

PO Box 709 Edison, NJ 08818 Phone: 908-769-5555 Fax: 908-222-1864 Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com