



Gobble up these amazing BLACK FRIDAY DEALS!

**BLACK
FRIDAY
SALE**
ALL MONTH
LONG


NEW

MC-24R™ REFRIGERATED MICROCENTRIFUGE

- * Touch control and digital display
- * Powerful refrigeration, as low as -20°C
- * High speed - 16,706 xg with COMBI-rotor
- * Quiet operation
- * Supplied complete with COMBI-Rotor™
- Accepts 24x1.5/2.0ml tubes AND 2xPCR strips

C2417-R-C	Refrigerated Microcentrifuge	\$5695.00
R2020	Roto-Mini™ Tube Rotator with adapters for common tube types	FREE! \$441.50
C2417-24SC	Rotor for 24 spin columns (17,317 xg)	\$820.30

FREE ROTO-MINI WITH PURCHASE - \$441.50 VALUE!



QuadCount™ Automated Cell Counter

- * User friendly, touch-screen interface
- * Auto-focus with motorized XYZ stage
- * LED light source, precision optics and CMOS detector
- * Quick, normal and precise counting modes
- * Process 4 samples at once with QuadSlide™ disposable slides



Just load the slide
and press start

E7500	QuadCount Automated Cell Counter	\$3151.00	4151.30
E7500-S1	QuadSlide 4 chambered disposable slides, pk/50		\$167.40

AUTOMATE YOUR CELL COUNTING! SAVE \$1000!

**Call Us Today!
BLACK FRIDAY
DEALS END
12/1/2023**

TC 9639 Gradient Thermal Cycler

- * Intuitive touch screen control
- * Program wizard - quick, easy and foolproof
- * Multiformat block for 96x0.2ml tubes, PCR strips, 39x0.5ml tubes or a 96 well plate
- * Advanced programming features such as adjustable ramp rate and time/temperature increments/decrements

T5000-96	TC 9639 Gradient Thermal Cycler	\$3500.00	\$5978.60
----------	---------------------------------	-----------	----------------------

HUGE DISCOUNT! SAVE OVER 40%!





Black Friday BOGO Deals

BOGO

BenchMixer™ V2

- * Reversible top: use "cup" or "flat/dimpled" surface
- * Powerful motor for instant vortexing
- * Q-Drive™ dynamic balancing system minimizes vibration
- * Variable speed with continuous or touch operation



BOGO - Put one on every bench!

BV1003 BenchMixer V2 with reversible top \$328.80/2

BUY ONE AT \$328.80 - GET ONE FREE!



BEADBUG MICROTUBE HOMOGENIZER

- * Powerful mixing action for lysis, grinding and homogenization
- * Eliminates cross contamination
- * Personal sized - 3 tube capacity
- * Cold room compatible
- * Wide variety of grinding media available



D1030 BeadBug 3, 3x2ml tube capacity \$980.80/2

D1032-30 TriplePure Zirconium Beads, 3mm, pk 50 tubes \$178.40

Additional bead materials/sizes available - see special offer below!

Perfect for small labs and student labs!

BUY ONE AT \$980.80 - GET ONE FREE!

Benchmark Homogenizer Bead Tubes

- * Convenient, prefilled, 2ml bead tubes
- * Compatible with most brands of bead beaters
- * Variety of bead sizes and materials to meet a range of sample needs
- * TriplePure™ line certified nuclease, protease and DNA free



D1032-30	TriplePure Bead Tubes, Prefilled, 3.0mm zirconium, pk 50 tubes	\$178.40
D1032-15	TriplePure Bead Tubes, Prefilled, 1.5mm zirconium, pk 50 tubes	\$178.40
D1032-10	TriplePure Bead Tubes, Prefilled, 1.0mm zirconium, pk 50 tubes	\$178.40
D1032-05	TriplePure Bead Tubes, Prefilled, 0.5mm zirconium, pk 50 tubes	\$178.40
D1032-01	TriplePure Bead Tubes, Prefilled, 0.1mm zirconium, pk 50 tubes	\$178.40
D1032-SK	TriplePure Starter Kit, 10 each of 0.1, 0.5, 1.0, 1.5 and 3.0mm	\$178.40
D1031-10	Prefilled Bead Tubes, 1.0mm glass, pk of 50 tubes	\$127.80
D1031-05	Prefilled Bead Tubes, 0.5mm glass, pk of 50 tubes	\$127.80
D1031-01	Prefilled Bead Tubes, 0.1mm glass, pk of 50 tubes	\$127.80
D1033-28	Prefilled Bead Tubes, 2.8mm steel, pk of 50 tubes	\$178.40
D1033-30G	Prefilled Bead Tubes, Garnet Shards, 1x6mm zirconium, pk/50	\$178.40

BUY 2 PACKS, GET ONE FREE - MIX OR MATCH!




BOGO

More BOGOs and Deals

Mini BlotBoy™

- * Tri-directional motion for thorough, gentle mixing
- * Optimized speed and tilt for blotting, washing and staining
- * 10.5 x 7.5" platform with rubber mat
- * Incubator and cold room safe

B3D1008	Mini BlotBoy 3D Rocker	\$485.30/2
B3D-STACK	Stacking platform, 3" clearance	\$135.80
B3D-DIMPLED	Dimpled mat for use with tubes	\$41.50



BUY ONE FOR \$485.30, GET ONE FREE!

myFUGE 12 Plus™


- * COMBI-Rotor™ for tubes and strips
- * Digital speed control - 7000rpm/2739 xg
- * Twice the capacity of traditional mini centrifuges
- * 99 minute timer
- * Near silent operation



Control Panel

C2212	myFuge 12 Plus Mini Centrifuge w/COMBI-Rotor	\$649.00/2
C1008-A2-6	Adapters for 0.2ml tubes, pk of 6	\$26.20
C1008-A5-6	Adapters for 0.5ml tubes, pk of 6	\$26.20

BUY ONE AT \$649, GET ONE FREE!



Top view of rotor

MULTI-THERM™

COOL - HEAT - SHAKE

- * Vortex mixer accommodates all common laboratory tubes and plates
- * Precision contoured wells for uniform thermal transfer
- * Exchangeable blocks for tubes 0.2ml to 50ml*
- * Two models: Heat/Cool or Heat only
- * Sequence function for linking or repeating programs (HC model only)

H5000-HC	MultiTherm Shaker with heating/cooling	\$3124.50
H5000-H	MultiTherm Shaker, heat only	\$2378.10
H5000-CMB	Combination block: 15x0.5ml & 20x1.5ml	\$340.60 FREE!

*Blocks sold separately. Additional blocks available
See www.benchmarkscientific.com for details

FREE COMBINATION BLOCK WITH PURCHASE - \$340 VALUE!



Check out these UNBELIEVABLE prices!



MHS | Series

Multi-Position Hotplate Stirrers

- * Heat and stir simultaneously
- * Digital control of time and speed
- * White, ceramic coated top surface
- * 3 positions
- * 1 Liter capacity at each position



IPS7000-3 3 Position Hotplate Stirrer ~~\$999.00~~ **\$699.00**
 *Larger models also available

SALE! SAVE 30%

SmartReader™ 96 Microplate Absorbance Reader

- * Stand alone operation - no computer needed
- * User-friendly, touchscreen programming
- * Integrated plate shaker
- * Supplied with 405, 450, 492 and 630nm filters

MR9600 SmartReader 96 Microplate Reader ~~\$6534.00~~ **\$5195.00**
P5208-200U ProPette™ 8 channel pipette, 20-200µl ~~\$399.00~~ **FREE**
MR9600-TP Thermal Printer for MR9600 ~~\$361.70~~
MR9600-PA Extra paper for printer, pk of 4 ~~\$34.50~~



SPECIAL PRICE AND FREE MULTICHANNEL PIPETTE - \$1738 VALUE!



MW9600 Microplate Washer


- * Precision needle position control
- * Auto-rinse to eliminate cross-contamination
- * Liquid level sensors in bottle
- * Customize and store up to 100 wash protocols

MW9600 Microplate Washer without manifolds ~~\$5009.30~~ **\$3999.00**

FREE 12 and 8 channel manifolds with purchase! FREE - \$709.60 value
MW9600-12CM 12 channel washer manifold ~~\$372.00~~ **FREE**
MW9600-8CM 8 channel washer manifold ~~\$337.60~~ **FREE**

SAVE \$1000 AND GET 2 FREE MAINFOLDS!

**AMAZING
BLACK
FRIDAY
DEAL!**



Free offers require purchase of designated items at listed price. Offers expire 12/01/2023. Not valid on previous purchases. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code TURKEY23.

Benchmark Scientific Phone 908-769-5555 Email info@benchmarkscientific.com