

Choose one
FREE

With purchase of item below*

(Free with items: C2000, B5200, H5100-HCT, T5 Series, B2410, MR9600/9610)

Your choice of:

*Promo Code BP1.24



Value: \$347

Includes microtube
& PCR rotors

MyFuge™ Mini

Item: C1008-G



Value: \$298

Speed Range:
200-3200 RPM

BenchMixer™

Item: BV1000



Value: \$100

Amazon Gift Card

Gift Card

Item: AMZ-100

PLATEFUGE™

MICROPLATE MICROCENTRIFUGE



- Quickly spins down droplets in PCR plates
- Unique rotor design prevents sample spillage
- 50% smaller than traditional centrifuges
- Accepts all popular PCR plates (skirted, semi & non) & most microtiter plates




Before/After Centrifugation During Centrifugation

C2000 PlateFuge™ Microplate Microcentrifuge **\$894**

SUREAIR™

PCR WORKSTATION WITH UV AND HEPA FILTRATION

- HEPA filter and UV light reduce cross-contamination during PCR preparation
- Compact, 55 x 40 cm footprint
- Cable port allows easy use of centrifuges, vortexers and other instruments
- Includes pipette holder and shelf for storage of pipettes and tips
- Provides clean air (to ISO 5 standards)



B5200 SureAir™ PCR Workstation **\$2,279**
Accessories:
B5200-FIL SureAir™ Replacement Filter **\$307**
B1450-IC UV-C Indicator Cards, 25 ct. **\$108**
B1450-ISS UV-C Indicator Stickers, 40 ct **\$96**

TC 9639

Thermal Cycler



- Intuitive touch screen control
- Gradient Capabilities
- Multi-format block accepts:
 - 0.2ml tubes or strips (96)
 - 0.5ml tubes (39)
 - 96 well plate (1)

Touch screen cycler
0.2ml tubes or strips (32)

T5000-96 TC 9639 Thermal Cycler, multi-format block **\$5,995**
T5000-384 TC 9639 Thermal Cycler, 384 well block **\$6,979**
T5005-3205 TC 32 Mini Thermal Cycler **\$3,775**

TC-32

Mini Thermal Cycler



BEAD BATH DUO™

THE ONLY BATHS DESIGNED FOR BEADS OR WATER



- Two separate modes for BEADS or WATER
- Quik-Cal™, simple in-lab calibration function
- Offers all of the advantages of a "Bead Bath"
 - Easy (accessory free) securing of tubes and bottles
 - Allows angled incubations
 - Runs constant (no wait for heat-up)



Choose bead or water mode

B2410 BeadBath Duo™ 10L **\$1,424**

MULTI-THERM™ Touch

COOL - HEAT - SHAKE



- Large color touch screen
- Multiple block options for plates & tubes up to 50ml
- Broad temperature & speed range
- Advanced programming features

Speed range: 200-3,000 rpm (block dependent)
Shaker Orbit: 3mm
Temp. range: Ambient -20° to 105°C



Heated lid prevents
condensation & evaporation

H5100-HCT MultiTherm Touch, Heating/Cooling Vortexer (blocks sold separately) **\$4,189**
H5100-HL Optional heated lid **\$803**
Available blocks*
H5100-05 24 x 0.5ml **\$374**
H5100-MP Block, microplate **\$374**
H5100-384 Block, 384 well plate **\$641**
*Additional blocks available, see website for details

SMARTREADER™ 96

Your go-to instrument for
ELISA and absorbance assays



- Intuitive, 7 inch, color touch screen control panel
- Filter based, wavelength range from 340-750nm (405, 450, 492, 630nm filter included)
- Operates as a stand alone system, with USB flash drive for data transfer
- Integrated, multi-speed plate shaker

Also Available:
MW9600
Microplate Washer

MR9600 SmartReader™ 96 Microplate Absorbance Reader **\$6,730**
MR9600-T SmartReader™ 96 Microplate Absorbance Reader with Incubation **\$7,164**
MR9600-TP SmartReader™ Thermal Printer **\$368**

BenchPress

Valid Through 07.31.24.
Pricing, availability, and specifications subject to change without notice.

Vol. 1.1 2024
Promo: BP1.24

Promos and Lab Favorites

Shakers Incubators Centrifuges Balances Rotators Vortexers DryBaths Rockers and MORE!



MC-24™ Touch

High Speed Microcentrifuge

- COMBI-Rotor™ accepts tubes and strips
- Full Color Touch Screen
- mySpin™ 1 button program recall
- High Speed, 16,800xg
- Rotor options for cryovials and 5ml tubes

Max. Speed: 13,500rpm / 16,800xg
Max Capacity: 24 x 1.5/2.0ml
2 x PCR Strips (16 x 0.2ml)



High Speed



Touch Control



Optional 8-place rotor for
5ml tubes and 2ml cryovials

C2417 MC-24™ Microcentrifuge with 24 place COMBI-Rotor **\$2,666**
C2417-R5 Optional 8 x 5ml rotor **\$1,128**
C2417-R5-ADP2 2ml cryovial adapter for C2417-R5 rotor, 6/pk **\$140**




C2212 myFuge™ 12 Plus with COMBI-Rotor **\$649**

myFUGE™ 12 Plus

Digital Mini Centrifuge

- COMBI-Rotor™ for tubes and strips
- Digital control up to 7,000 rpm / 2,739 xg!



C2417-R MC-24R™ Refrigerated Microcentrifuge with COMBI-Rotor **\$5,695**

MC-24R™

REFRIGERATED MICROCENTRIFUGE

- COMBI-Rotor™ accepts tubes and strips
- Powerful refrigeration, as low as -20°C
- Aerosol tight rotor options
- Status indicator light
- High Speed, 16,800 xg



High Speed Touch Control Powerful Refrigeration


QuadCount™

Automated Cell Counter

- Fast, accurate and reproducible counting & live dead determinations
- Intuitive software with auto focus and touch screen interface
- Cost-effective QuadSlides™ disposable slides with four individual chambers



E7500 QuadCount™ Automated Cell Counter **\$4,275**
E7500-S1 QuadSlide™ 4 chamber slides, pack of 50 **\$172**



Live/Dead cell (Viability)





BEADBLASTER™ 24 & 24R

HOMOGENIZER & REFRIGERATED HOMOGENIZER

- Completely lyse, grind & homogenize a variety of sample types
- Available with **REFRIGERATION** to protect delicate samples
 - True refrigeration - no dry ice or LN2 needed
 - Space saving design
- Process up to 24 samples simultaneously
- Wide variety of grinding media in prefilled 2ml & 5ml tubes

Also Available:



D2400-R	BeadBlaster 24R with refrigeration	\$15,098
D2400	BeadBlaster 24 Homogenizer	\$9,460
D2400-R5V2	Optional tube holder for 12X 5ml tubes	\$463
D1032-05	Prefilled 2ml tubes, 0.5mm zirconium, 50/pk.	\$182
D1032-15	Prefilled 2ml tubes, 1.5mm zirconium, 50/pk.	\$182
D1032-30	Prefilled 2ml tubes, 3mm zirconium, 50/pk.	\$182
D1033-28	Prefilled 2ml tubes, 2.8mm steel, 50/pk.	\$182
D1034-30	Prefilled 5ml tubes, 3mm zirconium, 50/pk.	\$244
D1034-28	Prefilled 5ml tubes, 2.8mm steel, 50/pk.	\$244

Additional grinding media is available in prefilled tubes and loose bulk. Ask for details.

Special Offer
\$13,895



Capacity: 24 x 2ml tubes
Speed range: 4.00 - 7.00 m/sec
D2400-R Only:
Cooling: Compressor
Temp. Range: -10°C to ambient



DryCHILL™

Ice-free Cooling Blocks

- Ice-free, stable cooling for sample tubes
- Unique design reduces condensation
- Maintains 4°C for up to 200 minutes
- 4 sizes for various tubes up to 50ml
- Visual temperature indicator



DC9602	DryChill™ 96 x 0.2ml/PCR Plate block	\$406
DC4020	DryChill™ 40 x 1.5/2.0ml block	\$406
DC1215	DryChill™ 12 x 15ml block	\$458
DC0650	DryChill™ 6 x 50ml block	\$458

INCU-SHAKER™ 20L/IR

STACKABLE SHAKING INCUBATOR

- Stainless steel chamber provides precise temperature uniformity
- Proven, maintenance-free Incu-Shaker™ motor drive system
- Intuitive touch-screen controls
- Included pull-out adjustable shelf for static incubations



MAGic-Clamp™ accessories for instant flask exchange

H3020	Incu-Shaker 20L with non-slip rubber mat	\$9,995
H3022	Incu-Shaker 20LR (refrigerated) with non-slip rubber mat	\$12,295

Stackable design for high throughput culturing



ROTO-THERM™

- Temperature control up to 60°C
- Includes tube holders for 0.5 to 50ml
- Easy, in-lab temperature calibration
- Available with variable speed & mix modes
 - Oscillating
 - Pausing
 - Rocking



Included Holders

- 11 x 1.5-2ml (2/pk.)
- 5 x 5-7ml (2/pk.)
- 4 x 10-15ml (2/pk.)
- 6 x 50ml* (1ea.)

*Also accepts 24 x 0.5ml

Temp Range: Ambient +5° to 60°C
Speed Range: 5 to 70rpm (Plus model)
Max. Capacity: 6 x 50ml, 16 x 15ml, 24 x 0.5ml, 44 x 1.5/2.0ml, 96 x 0.2ml

H2024	RotorTherm Plus w/included tube holders (variable speed)	\$889
R2020-1520	Optional tube holder for 96 x 0.2ml (8 x PCR strips)	\$88

PULSE Series

Ultrasonic Homogenizer

- Ideal for cell disruption, DNA shearing and ChIP assays
- Programmable with touch screen control and pulsed mode operation
- For volumes 0.1 to 150ml
- Supplied complete with soundproof enclosure, 6mm horn, temperature probe and integrated stand



DP0150	Pulse 150™ Ultrasonic Homogenizer with 6mm horn	\$3,531
DP0650	Pulse 650 Ultrasonic Homogenizer with 6mm horn	\$5,775
DP0150-2	Horn, 2mm diameter, for 0.1-5ml	\$325
DP0150-3	Horn, 3mm diameter, for 3-10ml	\$325
DP0150-8	Horn, 8mm diameter, for 25-150ml	\$325
DP0150-10	Horn, 10 diameter, 30-300ml samples, Pulse 650 only	\$325
DP0150-12	Horn, 12mm diameter, 50-500ml, Pulse 650 only	\$325

HERMLE Benchmark



5-Year Warranty!

Z446-K Refrigerated Universal Centrifuge

- Max. swing-out capacity: 28 x 50ml, 68 x 15ml, and 20 x microplates
- Temp Range -20°C to +40°C
- Bundle options for tubes or plates!
- Truly universal, 14 rotor options

Tissue Culture Bundle



Microplate Bundle

Z446-K-BND	Z446-K Refrigerated Tissue Culture Bundle	Reg. \$16,885	Bundle \$14,500
Z446-K-MPBND	Z446-K Refrigerated Microplate Bundle	Reg. \$17,055	Bundle \$14,100

INCU-SHAKER™

SHAKING INCUBATORS

- **MAGic™** Clamp accessories for instant-magnetic clamp exchange
- Available with CO₂ control or cooling
- **NEW!** Adjustable clamp for 125, 250 and 500ml flasks




Limited Time Offer*
FREE!
1 ea. **MAGic Clamp™** Universal platform for flasks & tube racks
PATENTED
Item: H1010-MR
H1000-MR

*Promo Code BP1.24

INCU-SHAKER Mini

Speed Range:	30 to 300rpm
Temp Range:	Amb. +5° to 70°C
Platform Dims:	9.5 x 11.5 in.
Universal Platform:	H1000-MR
Flask Capacity:	2 x 1L (H1000-MR-1000) 5 x 500ml (H1000-MR-500) 8 x 250ml (H1000-MR-250) 12 x 125ml (H1000-MR-125)

\$3,598



Item: H1001-M

INCU-SHAKER 10L & 10LR

Speed Range:	30 to 300rpm
Temp Range:	Amb. +5° to 60°C
CO ₂ Range:	Amb. -15° to 60°C
Platform Dims:	17.8 x 17.8 in.
Universal Platform:	H1010-MR
Flask Capacity:	5 x 2L (H1000-MR-2000) 9 x 1L (H1000-MR-1000) 16 x 500ml (H1000-MR-500) 25 x 250ml (H1000-MR-250) 36 x 125ml (H1000-MR-125)

**\$7,495/
\$10,596**



Item: H2010 & H2012 (Cooling)

INCU-SHAKER Mini CO₂

Speed Range:	30 to 300rpm
Temp Range:	Amb. +5° to 60°C
CO ₂ Range:	0 to 20%
Platform Dims:	9.5 x 11.5 in.
Universal Platform:	H1000-MR
Flask Capacity:	4 x 1L (H1000-MR-1000) 5 x 500ml (H1000-MR-500) 8 x 250ml (H1000-MR-250) 12 x 125ml (H1000-MR-125)

\$11,010



Item: H3501