


On **TARGET** 2024 Promotions

Vol. 1-2024
special offers valid
through Oct. 31, 2024



isoPURE™ Mini

Automated Purification System

- Compatible with all popular magnetic bead kits for DNA, RNA purification
- Process up to 16 samples
- Compact, 8x10 in. footprint
- Programmable by PC or Android app



Capacity:	16 samples
Plate Positions:	1 Plate Position
Lysis Temp. Range:	Ambient +5° to 120°C
Elution Temp. Range:	Ambient +5° to 120°C
Mixing Speed:	Variable (10 mixing levels)

AP1016 IsoPure™ Mini, Automated Purification System **\$10,447**

Also Available
isoPURE™ 96
Automated Purification System
\$47,250

Item:	AP1096
Capacity:	96 samples
Plate Positions:	8 Position carousel
Lysis Temp. Range:	Ambient +5° to 120°C
Elution Temp. Range:	Ambient +5° to 120°C
Mixing Speed:	Variable (10 mixing levels)



Special Offer

\$337
value



Purchase AP1016 and receive a free starter kit, including:

- One case AP1016-MC Tip combs
- One case AP1096-DWP Deep well plates



Smart-Q™ 100/800

Quantitative Fluorometer For quantification of DNA, RNA & Protein

- Sensitive quantitation down to 0.1pg/μl ds DNA
- Large color touch screen for programming and operation
- Two models for single or 8-sample processing
- Built-in reagent calculators for sample preparation
- Compact and portable (11 x 6 inch footprint)

CALCULATE RESULTS
IN SECONDS

SINGLE OR
MULTICHANNEL

Display:	7.0-inch, color touch screen
Dynamic Range:	5 orders of magnitude
Light Source:	Blue and Red LED
Consumable type:	0.5mL assay tubes


Choose 1 FREE Quantification Kit with purchase of F1001 OR F1008

F1000-HS-100	High Sensitivity dsDNA Quant. Kit
F1000-BR-100	Broad Range dsDNA Quant. Kit

Value: \$131

Use promo code: **OT1.24**

F1001	Smart-Q™ 100 Fluorometer for single samples	\$3,900
F1008	Smart-Q™ 800 Fluorometer for 8-samples	\$6,650



SMART

SAFE Nucleic Acid Stain

GLOW

- Replaces hazardous EtBr
- Excitation by UV or blue light

Use promo code: **OT1.24**


E4500-LD	SmartGlow™ Loading Dye, 1.0ml	\$124
E4500-PS	SmartGlow™ Safe Green Pre-Stain, 1.0ml	\$124

FREE

Mini Blue Light
Transilluminator
with purchase of
2 tubes of
SmartGlow™ stain

Value:

\$491



E4100

Accuris Analytical Balances

- Readability of 0.0001g
- Available with Quick-Cal™ internal calibration feature
- Multiple mode options (standard weighing, count mode, & percent)

Analytical Balances

W3100-120	120g x 0.0001g	\$2,143
W3100A-120	120g x 0.0001g w/ QuickCal™	\$2,274
W3100-210	210g x 0.0001g	\$2,488
W3100A-210*	210g x 0.0001g w/ QuickCal™	\$2,637

Compact Balances

W3300-120	120g x 0.01g	\$327
W3300-500	500g x 0.01g	\$381
W3300-1200	1200g x 0.01g	\$424

***FREE!**
With purchase of
W3100A-210

***FREE!**
Compact
Balance

With purchase of
W3100A-210

Accuris Compact Balances

- Readability of 0.1g to 0.01g
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer

Use promo code: OT1.24



SmartDrop™ Nano Spectrophotometers

- Quick and accurate results from a sample droplet (0.5µl - 2.0µl)
- 7.0-inch color touch screen for operation (no PC software required)
- Multifunctional software for Nucleic Acids, Protein A280, Colorimetry, UV-Vis, & OD600 measurements
- Built-in printer and cuvette port for OD600 readings
- Three models available (SmartDrop™ L, X & XF)

NS1000	SmartDrop™ L	260,280,600nm	\$6,090
NS1010	SmartDrop™ X	200-800nm	\$11,539
NS1020	SmartDrop™ XF	200-800nm with fluorometer	\$13,639


BUNDLE OFFERS: Special price, includes FREE pipette

NS1010-BNDL SmartDrop™ X Bundle \$8,000
NS1020-BNDL SmartDrop™ XF Bundle \$9,400

**SmartDrop™ XF Model
Includes
Fluorometer!
(for DNA quantification)**

***FREE!**
NextPette™
Pipette
with purchase
of bundle

**P7700-10
0.5-10µl
\$276 value**
Use promo code: OT1.24



QuadCount™ Automated Cell Counter

- Reduce counting errors, save time and improve efficiency
- Intuitive software and touch-screen interface
- Utilizes cost-effective, disposable, 4-chamber QuadSlides™
- Quick, normal and precise count modes (up to 3.6µl count volume)

E7500 QuadCount™ \$4,275
E7500-S1 QuadSlide™, 50/pk. \$172

Sale Price: \$4,075

Use promo code: OT1.24



Replaces manual cell counting...ask for a demo today!


**SAVE
\$200**



NEW PRODUCT

Accuris UV Transilluminator

- Large, 21x25cm viewing surface
- Powerful 14 watt bulbs
- Variable intensity
- 302nm midrange UV wavelength
- Optimal for EtBr, SmartGlow™, and SYBR stains
- Hinged, UV-blocking cover



E3200 UV Transilluminator, Variable Intensity **\$1,250**

**Intro Price
\$1,250**