

hybex™

Media Storage Bottles



**Borosilicate
Glass**

- Borosilicate 3.3 hybex glass
- Autoclavable caps (GL45) & drip rings
- Permanent enamel graduation markings
- Large area for bottle marking/identification
- Autoclavable, reuseable & breakage resistant
- Five cap colors

Technical Data

Material:	hybex™ borosilicate glass (3.3)
Cap Thread Size:	GL45*
Temp. Range:	-70°C to 500°C
Sterilization:	Autoclave or dry heat
Cap Material:	Polypropylene (140°C)
Cap Colors:	Blue, green, purple yellow, white or assorted

*50ml bottles are only supplied with blue GL32 sealing caps

Bottle Dimensions

B3000-50	46 x 87 mm
B3000-100	56 x 100 mm
B3000-250	70 x 138 mm
B3000-500	86 x 176 mm
B3000-1000	101 x 225 mm
B3000-2000	136 x 260 mm
B3000-5L	330 x 187 mm
B3000-10L	410 x 235 mm

Also Available

**hybex™
Erlenmeyer Flasks**



Ordering Information

B3000-50	Media storage bottle, 50ml with standard (GL32) blue cap, 10/pk.
B3000-100*	Media storage bottle, 100ml with standard (GL45) cap, 10/pk.
B3000-250*	Media storage bottle, 250ml with standard (GL45) cap, 10/pk.
B3000-500*	Media storage bottle, 500ml with standard (GL45) cap, 10/pk.
B3000-1000*	Media storage bottle, 1L with standard (GL45) cap, 10/pk.
B3000-2000*	Media storage bottle, 2L with standard (GL45) blue cap, 5/pk.
B3000-5L	Media storage bottle, 5L with standard (GL45) blue cap, 1/ea.
B3000-10L	Media storage bottle, 10L with standard (GL45) blue cap, 1/ea.
B3000-SP*	Media storage bottle, Starter Pack with standard (GL45) cap (Starter Pack of ten bottles includes 2x100ml, 3x250ml, 3x500ml & 2x1000ml)

B3000-HTC	High temperature red caps, up to 180°C (GL45), 10/pk
B3000-CAP*	Replacement cap, standard blue (GL45), 10/pk.
B3000-RIN	Replacement sealing ring, (GL45), 10/pk.
B3000-CAP2	Replacement cap, standard blue (GL32), 10/pk.
B3000-RIN2	Replacement sealing ring, (GL32), 10/pk.

*Please add (-B) for blue caps, (-G) for green caps, (-P) for purple caps, (-Y) for yellow caps, (-W) for white caps or (-AS) for assorted caps (2 of ea. color)

Benchmark
Scientific



Benchmark Scientific Inc.
PO Box 709 Edison, NJ 08818
Phone: 908-769-5555 Fax: 908-222-1864
Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com