

Celebrate the Luck O'the Irish St. Patrick's Day DEALS!



St. Patrick's Day Sale

March 10 - April 11
Offer code
SHAMROCK25



myFUGETM Mini

- * Complete with microtube rotor, strip tube rotor and adapters for 0.5ml tubes
- * Click-n-LokTM rotor exchange system - no tools!
- * Improved braking, decelerates in just 1 second!
- * Big capacity in small space:
 - ✓ 8 x 1.5/2.0ml or 0.5ml tubes
 - ✓ 4 x PCR strips
 - ✓ 32 x 0.2ml tubes

C1008-G myFuge Mini, 6000rpm/2000xg, green lid \$346.90/2

BUY ONE AT \$346.90 - GET ONE FREE!



On-board storage in base for PCR rotor and adapters

ACCURIS UV Transilluminator



- * High powered, 14 watt bulbs for superior sensitivity
- * Large viewing surface - 21x26cm
- * Hinged, UV-blocking cover
- * Continuously adjustable UV intensity
- * 302nm - ideal for Ethidium Bromide, SmartGlowTM, SYBR[®] Green and other common stains

E3200 AccurisTM UV Transilluminator, 21x26cm ~~\$1250.00~~
SALE PRICE \$995.00

SAVE SOME GREEN - \$255 OFF!



myTEMP^{Mini} CO₂TM

- * Personal sized - 13x14.5" footprint (20L)
- * Digital control of temperature and CO₂
- * Dual beam, IR sensor
- * Accommodates bottles and flasks up to 2L
- * Supplied with 2 stainless steel shelves

H2300-HC2

myTemp Mini CO₂ Digital Incubator

Sale Price
~~\$3695.00~~ \$4586.00



KEEP YOUR CELLS HAPPY AND SAVE!

FREE STUFF AND AMAZING DEALS!



**FREE
Agarose!**



- * Excellent replacement for Ethidium Bromide
- * Excitation by UV or blue light
- * Non-mutagenic, environmentally safe
- * Available as a pre stain (PS) or loading dye (LD)

SMART
SAFE Nucleic Acid Stain
GLOW

Turn your samples **GREEN!**
Free 100g Agarose with purchase of TWO SmartGlow™!!

E4500-PS	SmartGlow Pre stain, 1ml (20,000X)	\$124.95
E4500-LD	SmartGlow Loading Dye with stain, 1ml (6X)	\$124.95
A1701	Agarose LE, General Purpose, 100g	FREE! \$139.87



FREE! 100G AGAROSE WITH SMARTGLOW!



VORNADO™
MINIATURE VORTEX MIXER

- * Powerful vortexing for tubes up to 50ml
- * "Instant-touch" pressure activated operation
- * Mini sized - fits in the palm of your hand!
- * Shape of mixing cup protects motor from spills



Call Us Today!
THE LUCK RUNS
OUT 4/11/2025



BV101-G Vornado Mini Vortexer, green head \$229.09/2

BUY ONE AT \$229.09, GET ONE FREE - \$229.09 VALUE

LIMITED TIME ONLY
**BUY 2 GET 1
FREE**

Benchmark Homogenizer Bead Tubes

- * Convenient, prefilled, 2ml bead tubes
- * Compatible with most brands of bead beaters
- * Variety of bead sizes and materials to meet a range of sample needs
- * TriplePure™ line certified nuclease, protease and DNA free

D1032-30	TriplePure Bead Tubes, Prefilled, 3.0mm zirconium, pk 50 tubes	\$184.69
D1032-15	TriplePure Bead Tubes, Prefilled, 1.5mm zirconium, pk 50 tubes	\$184.69
D1032-10	TriplePure Bead Tubes, Prefilled, 1.0mm zirconium, pk 50 tubes	\$184.69
D1032-05	TriplePure Bead Tubes, Prefilled, 0.5mm zirconium, pk 50 tubes	\$184.69
D1032-01	TriplePure Bead Tubes, Prefilled, 0.1mm zirconium, pk 50 tubes	\$184.69
D1032-SK	TriplePure Starter Kit, 10 each of 0.1, 0.5, 1.0, 1.5 and 3.0mm	\$184.69
D1031-10	Prefilled Bead Tubes, 1.0mm glass, pk of 50 tubes	\$132.30
D1031-05	Prefilled Bead Tubes, 0.5mm glass, pk of 50 tubes	\$132.30
D1031-01	Prefilled Bead Tubes, 0.1mm glass, pk of 50 tubes	\$132.30
D1033-28	Prefilled Bead Tubes, 2.8mm steel, pk of 50 tubes	\$184.69
D1033-30G	Prefilled Bead Tubes, Garnet Shards, 1x6mm zirconium, pk/50	\$184.69



BUY 2 PACKS, GET ONE FREE - MIX OR MATCH!

MORE FREE STUFF! SAVE SOME GREEN!



MORTEXER™

MORE Power, MORE Flexibility... MORE Vortexer!

- * Unique Multi-Head™
-hands free microtube vortexing
- * Instant vortexing of tubes up to 50ml
- * Q-Drive™ dynamic balancing system minimizes vibration
- * Continuous or touch operation

BV1005 Mortexer Vortex Mixer w/Multi-Head \$341.70/2



Vortexing a 50ml tube in the Mortexer™



BUY ONE AT \$341.70 - GET ONE FREE!

MY GEL™ InstaView



- * View migration as it happens
- * Blue light illumination will not damage DNA
- * Take high resolution images without removing gel from the apparatus
- * Complete electrophoresis solution:
 - Gel box and casting trays
 - Direct connect power supply
 - Blue LED transilluminator
 - Imaging enclosure for taking gel pictures with a cellphone or tablet



E1201 myGel InstaView Complete Electrophoresis System \$1,016.73
 E4001-VG1 SmartBlue™ Viewing Glasses (for green dyes/blue light) \$53.48
 A1701 Agarose LE, General Purpose Agarose, 100g bottle ~~\$139.87~~ **2 FREE!**



TWO FREE BOTTLES OF AGAROSE - \$279.74 VALUE



Taq MasterMix and SmartCheck™ DNA Ladders

- * Ideal for routine PCR, genotyping, colony PCR and fast PCR
- * Improved template affinity for higher enzyme activity and yield
- * Hot Start version ideal for low copy number templates
- * Optimized proprietary buffer
- * 2X mastermix with inert red loading dye



- * Three sizes meet most applications
- * Ready-to-use with buffer and tracking dye
- * 100 lanes per tube (5µl per lane)
- * Higher intensity reference bands
- * Ultra-pure production allows for room temperature shipping

PR1001-R-1000	Taq 2X MasterMix, w/red dye, 1,000x50µl rxns	\$383.36	PR4005-500	50bp DNA Ladder, 500µl	\$189.53
PR1001-HSR-200	Hot Start Taq 2X MasterMix, w/red dye, 200x50µl rxns	\$196.36	PR4010-500	100bp DNA Ladder, 500µl	\$189.53
A1701	Agarose LE, General Purpose, 100g	\$139.87 FREE!	PR4100-500	1kb bp DNA Ladder, 500µl	\$189.53
			A1701	Agarose LE, General Purpose, 100g	\$139.87 FREE!



FREE AGAROSE WITH PURCHASE! - \$139.87 VALUE!

CATCH THESE DEALS BEFORE THEY'RE GONE!

myBLOCK™ Digital Dry Baths

- * Hinged cover improves uniformity and conserves energy
- * Optional external temperature probe for control of actual sample temperature
- * Unique Quick-Flip™ universal micro block accepts: 0.5ml, 1.5ml, 2ml, and 0.2ml tubes or 96 well plate*

Purchase a myBlock I and get a **FREE** Quick-Flip Block
Purchase a myBlock II and get **2 FREE** Quick-Flip Blocks

BSH5001	myBlock I, single chamber, w/o blocks	\$554.82
BSH5002	myBlock II, dual chamber, w/o blocks	\$690.75
BSWCMB	Quick-Flip universal micro-block FREE!	\$203.92
BSH-TP1	Optional external temperature probe	\$239.54

Additional blocks available. Visit website for details.
*96 well plate requires myBlock II and 2 Quick-Flip blocks



FREE QUICK-FLIP BLOCK WITH PURCHASE - UP TO \$407 VALUE

SureTemp™ Ovens

- * Temperature range to 200°C
- * Stainless steel chamber
- * Fiberglass insulation
- * Large viewing window
- * In-lab calibration feature
- * Three sizes: 30L, 50L and 150L

		Sale Price	Sale Price
H2606-30	Benchmark SureTemp Oven, 30L	\$1295.00	\$1095.00
H2606-50	Benchmark SureTemp Oven, 50L	\$1595.00	\$1295.00
H2606-150	Benchmark SureTemp Oven, 150L	\$2895.00	\$2495.00

SAVE \$100 ON A SURETEMP OVEN!



TC 9639 Gradient Thermal Cycler

- * Intuitive, touch screen control
- * Program Wizard - quick, easy and foolproof
- * Unique multiformat block - accepts 0.2ml tubes or strips, 96 well plate or 0.5ml tubes
- * Fast and reproducible results

		Sale Price	Sale Price
T5000-96	TC 9639 Gradient Thermal Cycler	\$4,495.00	\$5,995.00
C2000	PlateFuge™ for 2 microplates	FREE!	\$894.54

🍀 Ask about our full line of PCR Reagents

SAVE \$1,500 AND GET A FREE PLATEFUGE!



Free offers require purchase of designated item at listed price. Offer expires 4/11/2025. Not valid on previous purchases. Special price applies only to item specified. Valid in the USA only. Use promo code SHAMROCK25.