

Roto-Bot Mini

Roiú-Bot

	cito-Bot Mini			Roto-Bot <sup>™</sup>	
Tube Holder Capacity & Ordering Information:					
Item Number:	R4040-1520	R4040-150	R4040-500	R4040-HZ15	R4040-HZ50
Description:	6 x 1.5ml	1 x 15ml	1 x 50ml	6 x 15ml	6 x 50ml
Tube Diameter:	10 to 13mm	16 to 17mm	28 to 30mm	16 to 17mm	28 to 30mm
Orientation:	Dual*	Vertical	Vertical	Horizontal	Horizontal
Holder Capacity (Mini/Standard):	10/16	10/22	6/12	1/2	1/2
Tube Capacity (Mini/Standard):					
Package Unit:	60/96 Pack of 2	10/22 Pack of 2	6/12 Pack of 2	6/12 1 ea.	6/12 1 ea.

The RotoBot<sup>™</sup> Rotating Mixers feature advanced programming options that far surpass that of the traditional laboratory rotator. With digital speed control and more than 10 operating modes, the RotoBot can be selected to perform a wide variate of rotation of these actions can be selected to perform a wide variate of rotation. providing the user with a nearly infinite level of programming.

The RotoBot and RotoBot Mini each inlcude an assortment of vertical tube holders. The Mini model is supplied with holders for 12 x 1.5ml tubes, 2 x 15ml & 2 x 50ml tubes; the standard model, 30 x 1.5/2.0ml, 8 x 15ml & 2 x 50ml tubes. For higher capacity applications, more of these holders can be purchased as accessories. In addition, horizontal holders are available for applications such as Western blotting and hybridization. The unique, sliding attachment allows for quick and easy changing of tube holders.



For rocking applications, a platform with a non-slip rubber mat is available. Once attached, simply select the "PP" mode to turn the rotator into a platform rocker.

F1	F2	F3	F4	F5	F6	F7	F8	F9	 C1 C2	PP*	Key:
$\bigcirc$	STOP	STOP STOP	* STOP	STOP	$\sum$	STOP STOP	Carop Store	Ç	Custom Mode	Platform Rocking Mode	Pause Mix

## Technical Data:

Speed Range: Capacity (R4040): Capacity (R4045): Timer: W x D x H (R4045) **Operating Temp. Range:** Weight: **Electrical:** 

1 to 65rpm 6 x 50ml, 10 x 15ml or 60 x 1.5ml 12 x 50ml, 22 x 15ml or 84 x 1.5ml 1 min. to 99 hours 59 min. / Cont. +2°C to +70°C 4 lbs/1.8 kg (Mini) / 5 lbs/2.27kg (St 100 to 240V, 50-60Hz, 3A

## **Ordering Information:**

	R4040*	<b>RotoBot Mini</b> with tube holder starter set Includes 6x1.5ml holders (2ea.) 1x15ml holder (2ea.) & 1 x 50ml holder (2ea)
cm	R4045*	<b>RotoBot with</b> with tube holder starter set Includes 6x1.5ml holders (5ea.) 1x15ml holder (8ea.) & 1x50ml holder (2ea)
cm td.)	R4040-RP R4045-RP2	Rocking Platform (w/ flat mat) for R4040, 6 x5.5 ir Rocking Platform (w/ flat mat) for R4040, 6 x10 in

Includes US plug. For EU plug, please add (-E)

Benchmar Scientific

**PO Box 709 Edison, NJ 08818** Phone: 908-769-5555 Fax: 732-313-7007 Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com