

**BIG SAVINGS**  
Available  
through  
**September 30<sup>TH</sup>**  
**2025**

# BACK TO SCHOOL 2025 PROMO

Kickstart the Semester with Lab Favorites for **\$2,025**



## Accuris™ TX Balance, 220g

Item: W3102A-220 ~~\$2,977~~ **\$2,025**

- **Electromagnetic Force Restoration (EFR) technology, fast response and stabilization**
- **Full color, touch screen control panel**
- **In-lab calibration function**

Capacity: **220 grams**  
Readability: **0.0001g**  
Pan size: **80mm**



Unique Combi-Rotor  
Accepts tubes and strips



Full Color Touchscreen

## MC-24™ Touch High Speed Microcentrifuge

Item: C2417 ~~\$2,693~~ **\$2,025**

- **Full color touch screen**
- **COMBI-Rotor™ accepts tubes and strips**
- **Unique air-flow system, keeps samples cool**
- **Up to 16,800 xg!**

Speed Range: **200 to 13,500 rpm (up to 16,800xg)**  
Capacity: **24 x 1.5/2.0ml**  
**2 x PCR Strips (16 x 0.2ml)**



## Incu-mixer™ MP Heated Microplate Vortexer, 4 position Mixer

Item: H6004 ~~\$2,961~~ **\$2,025**

- **Optimized for incubated mixing in microplates**
- **Heated lid prevents condensation**
- **Accepts plates up to 46mm high**
- **Independently control time, temp. and speed**

Temp. Range: **Ambient +5 to 70°C**  
Shaking Orbit: **2mm**  
Shaking Speed: **100-1,500 rpm**



## MultiTherm™ Thermal Shaker

Item: H5000-HC ~~\$3,363~~ **\$2,025**

- **Heats, cools and shakes**
- **Exchangeable blocks, for tube sizes 0.2 to 50 ml**
- **Precision contoured wells for uniform thermal transfer**

H5000-05	Block, 54 x 0.5ml	\$368
H5000-15	Block, 35 x 1.5ml	\$368
H5000-20	Block, 35 x 2.0ml	\$368
H5000-150	Block, 12 x 15ml	\$519
H5000-500	Block, 6 x 50ml	\$519
H5000-02	Block, 96 x 0.2ml or 1 x PCR plate	\$437



## SureAir™ PCR Workstation

Item: B5200 ~~\$2,383~~ **\$2,025**

- **UV-C light destroys stray DNA**
- **HEPA filter provides Class 5 cleanroom conditions**
- **Cable access port allows for powering small equipment in the chamber**

B1450-IS	UVC indicator sticker, 40/pk	\$102
B1450-IC	UVC indicator cards, 25/pk	\$115
B5200-FIL	Replacement chamber filter	\$322
B5200-PRE	Replacement pre-filter	\$52

Light Source: **8W UV-C Bulb**  
Int. Dimensions: **19.7 x 15.4 x 15.7 in.**



## BeadBug™ 6 Six Position Homogenizer

Item: D1036 ~~\$2,982~~ **\$2,025**

- **Extremely powerful mixing for lysis, grinding or homogenization**
- **Faster and more effective than tissue grinders**
- **Eliminates the cross contamination associated with handheld homogenizers**

D1032-10	Triple-Pure Zirconium, Ø 1.0mm, 2ml	\$184
D1032-15	Triple-Pure Zirconium, Ø 1.5mm, 2ml	\$184
D1032-30	Triple-Pure Zirconium, Ø 3.0mm, 2ml	\$184
D1033-28	Stainless steel, Ø 2.8mm, 2ml	\$184
D1034-28	Stainless Steel, Ø 2.8mm, 5ml	\$246

Speed Range: **2500 to 4350 rpm**  
Capacity: **6 x 2ml or 2 x 5ml**




Speed Range: **200 to 1500 rpm**  
Temp. Range: **Ambient -20 ° to 100°C**  
Orbit: **3mm**



# Also Available from Benchmark:

Valid Thru. 9.30.25  
Promo Code: BTS1

## myGel™ InstaView Mini Electrophoresis System



- View DNA migration and separation as it happens
- Blue light illumination – won't damage DNA
- Capture high resolution images without removing the gel from the tank
- Compatible with all smart phones for imaging

Item: E1201: ~~\$1,108~~ **\$995**

Output Voltage: 35V, 50V or 100V  
Gel Capacity: One 10.5x10cm gel, 0.5cm thick  
Timer: 0-99 minute or continuous

## Accuris™ UV Transilluminator

- 302nm midrange UV wavelength
- Scratch-resistant, UV filter glass
- Optimal for EtBr, SmartGlow™, and SYBR® Stain
- Hinged, UV-blocking safety cover



Light Source: 4 x 14-watt bulbs  
Intensity Control: Full variable intensity  
Viewing Surface: 21 x 26cm

Item: E3200: ~~\$1,250~~ **\$995**

## myTemp™ Mini CO<sub>2</sub> Digital Incubator



- Digital control of temperature and CO<sub>2</sub>
- Personal sized, 13 x 14.5" footprint (20L)
- Accepts bottles and flasks up to 2L
- Includes 2 stainless steel shelves
- Dual beam, infrared sensor

Item: H2300-HC2: ~~\$4,678~~ **\$3,995**

Temperature Range: Ambient -10°C to 55°C  
Temperature Accuracy: 0.25°C (at 37°C)  
CO<sub>2</sub> Range: 0 to 20%

## Orbi-Shaker™ CO<sub>2</sub> Remote Control Shaker

- Remote controlled, safe for use in CO<sub>2</sub> environments
- Horizontally circular 19mm orbit for aeration & mixing
- Instantly exchange flask clamps with MAGIC Clamp™ accessories (PATENTED)



Speed Range: 30 to 300 rpm  
Shaking Orbit: 19mm  
Operating Temp. Range: +0°C to +65°C

Item: BT4001: ~~\$3,879~~ **\$2,995**

## MC-24R™ Refrigerated Microcentrifuge



- Powerful refrigeration, as low as -20°C
- Touch operation
- Up to 16,800 xg!
- Included COMBI-Rotor™ accepts tubes and strips (lid included)

COMBI-rotor  
24 x 1.5/2.0ml  
AND  
2 x PCR strips



Item: C2417-R: ~~\$5,924~~ **\$4,995**

Max Speed: 13,500 rpm / 16,800xg  
Temperature Range: -20°C to +40°C  
Capacity: 24 x 1.5/2.0ml  
4 x PCR Strip (32 x 0.2ml)

## SmartWasher™ 96 Plate Washer


- Precision control of aspiration needle position
  - Accurate to 0.1mm
  - Minimal residual liquid, ≤1µl/well
- Auto rinse function eliminates cross contamination
- Up to 100 wash protocols can be customized and stored



Wash Heads: 8 channel or 12 channel (purchased separately)  
Wash Cycles: 1 to 99 cycles  
Program Storage Capacity: 100 programs

Item: MW9600: ~~\$5,472~~ **\$4,995**

## Incu-Shaker™ 10L Shaking Incubator



- TOUCH screen control of time, speed and temperature
- Exceptional temperature uniformity
- Horizontal circular 19mm orbit for aeration & mixing
- Instantly change flask clamps with MAGIC Clamp™ accessories
- Large platform 18 x 18", up to 5 x 2L

Item: H2010: ~~\$7,570~~ **\$5,995**

Temperature Range: Ambient +5° to 60°C  
Speed Range: 30 to 300 rpm  
Shaking Orbit: 19mm (3/4")

## SmartReader™ 96 Absorbance Plate Reader

- Intuitive, 7 inch color touch screen
- control panel
- Filter based system with wavelength range from 340 to 750nm
- Operates as a stand alone



Item: MR9600: ~~\$7,103~~ **\$5,995**

Light Source: Quartz-halogen lamp 6V/10W  
Readout Range: 0 to 4.000Abs  
Test Speed: Less than 6 sec.