

BLACK FRIDAY

SALE

Mini Mixer™



- Includes dimpled mat; ideal for mixing in blood & centrifuge tubes
- Tri-directional motion for thorough, uniform mixing
- Optional stacking platforms available
- Cold room and incubator safe

**B3D1020
C1008**

**MiniMixer™ 3D Nutating Shaker
MyFuge™ Mini Centrifuge**

\$473

~~\$347~~ FREE!

B3D-STACK-D

Stacking platform, small 10.5"x7.5" with dimpled mat

\$154

B3D-FLAT

Flat Mat, small 10.5" x 7.5"

\$47

B3D-DIMPLED

Dimpled Mat, small 10.5" x 7.5"

\$44

Specify lid color when ordering myFuge™ Mini - use -B for blue, -G for green, -P for purple, -C for clear and -R for red

**FREE MYFUGE MINI WITH PURCHASE!
\$347 VALUE**

ACCURIS COMPACT BALANCES



- User friendly interface
- Easy in-lab calibration
- 0.01g readability
- RS232 interface for data output

W3300-500

Accuris Compact Balance, 500g capacity

\$403

C1008

myFuge™ Mini Centrifuge

~~\$347~~ FREE!

Specify lid color when ordering myFuge™ Mini - use -B for blue, -G for green, -P for purple, -C for clear and -R for red

FREE MYFUGE™ MINI WITH PURCHASE!

myTEMP™ Mini Digital Incubators

- Digital temperature control
- Personal sized, 13 x 14.5" footprint (20 L)
- Heats and cools!
- Accepts bottles and flasks up to 2 L

H2200-HC myTemp™ Mini, with heating and cooling \$962/2



**BUY ONE AT
\$962, GET ONE
FREE!**

VORNADO™ MINIATURE VORTEX MIXER

- Powerful vortexing for tubes up to 50ml
- Inst-Touch™ pressure activated operation
- "Mini sized", less than 4 in. wide
- Available in 5 colors

BV101 Vornado™ Mini Vortex Mixer \$229/2
with gray cup head



**BUY ONE
AT \$229,
GET ONE
FREE**



Specify cup color when ordering- use -B for blue, -G for green, -P for purple, and -R for red


BENCHMARK SCIENTIFIC PHONE (908)769-5555 EMAIL: INFO@BENCHMARKSCIENTIFIC.COM

BLACK FRIDAY DEALS



myFUGE™ 12 Plus

- COMBI-Rotor™ for tubes and strips
- Digital speed control - 7000rpm/2739 xg
- Twice the capacity of traditional mini centrifuges
- 99 minute timer
- Near silent operation



DIGITAL DISPLAY




COMBI-ROTOR

C2212 myFuge™12 Plus Digital Mini Centrifuge **\$678/2**
BUY ONE AT \$678. GET ONE FREE



PLATEFUGE™ MicroPlate MicroCentrifuge

- Quickly spins down droplets in PCR plates
- Unique, patented swing-out rotor design prevents spillage
- Half the size of traditional centrifuges
- Accepts all popular types of PCR plates (skirted, semi, and non) and most microtiter plates



BEFORE/AFTER
CENTRIFUGATION



DURING
CENTRIFUGATION

C2000 PlateFuge™ MicroPlate MicroCentrifuge **\$895**

C1008 MyFuge™ Mini Centrifuge

\$347 FREE!

C2000-A02 Adapter for 0.2ml tubes and strips

\$67

C2000-DWMP Optional rotor for deepwell plates

\$236

Specify lid color when ordering myFuge™ Mini - use -B for blue, -G for green, -P for purple, -C for clear and -R for red

FREE MYFUGE MINI WITH PURCHASE!
\$347 VALUE



BEADBUG™ MicroTube Homogenizer

- Powerful mixing action for lysis, grinding and homogenization
- Eliminates cross contamination
- Personal sized - 3 tube capacity
- Cold room compatible
- Wide variety of grinding media available



**BUY
ONE AT
\$998,
GET ONE
FREE**

D1030 BeadBug™ 3, 3x2ml tube capacity

\$998/2

Prefilled Bead Tubes

D1032-30 TriplePure Zirconium Beads, 3mm, pk 50 tubes

\$185

D1032-15 TriplePure Zirconium Beads, 1.5mm, pk 50 tubes

\$185

D1032-10 TriplePure Zirconium Beads, 1.0mm, pk 50 tubes

\$185

D1032-05 TriplePure Zirconium Beads, 0.5mm, pk 50 tubes

\$185

D1032-SK TriplePure Starter Kit, 10 each 0.1, 0.5, 1.0, 1.5 and 3.0mm **\$185**

Other bead material/sizes available. Ask for details



BLACK FRIDAY DEALS

DIGITAL HEAT BLOCK



- Digital accuracy and convenience
- Assorted blocks for a wide variety of tube sizes
- Custom blocks available upon request
- Blocks sold separately



\$181 OFF, DON'T WAIT ANY LONGER!

BSH1001	Benchmark Digital Dry Bath, single position, without blocks	\$580 \$399
BSW1500	Block, 24 x 1.5ml	\$158
BSW1520	Block, 24 x 1.5/2.0ml	\$158
BSW05	Block, 24 x 0.5ml	\$158
BSW15	Block, 12 x 15ml (block height is 3.25in.)	\$158
BSW50	Block, 5 x 50ml (block height is 3.25in.)	\$158
BSW10	Block, 20 x 10mm (or 20 x 2.0ml)	\$158
BSW5MT	Block, 12 x 5ml centrifuge tubes	\$158
BSW1516	Block, 12 x 15-16mm (or 12 x 10ml)	\$158

*Additional blocks available




BSW1516



BSW1520

INCU-SHAKER™ 10L



- Large platform 18 x 18, up to 5 x 2 L
- TOUCH screen control over time, speed and temperature
- Exceptional temperature uniformity
- Horizontal circular 19mm orbit for aeration & mixing
- Instantly change flask clamps with MAGic Clamp™ accessories*
- Optional shelf for Petri dishes, microplates, etc.

OVER 30% OFF! LIMITED TIME ONLY!

H2010	Incu-Shaker™ 10L with non-slip rubber mat	\$7,570 \$4,995
H1010-MR	MAGic Clamp universal platform	\$439
H1000-MR-CMB	MAGic Clamp for 125ml, 250ml, or 500ml flask	\$56
H1000-MR-125	MAGic Clamp, 125ml	\$48
H1000-MR-250	MAGic Clamp, 250ml	\$48
H1000-MR-500	MAGic Clamp, 500ml	\$48
H1000-MR-1000	MAGic Clamp, 1000ml	\$73
H1000-MR-2000	MAGic Clamp, 2000ml	\$95

*Additional clamp sizes available. MAGic Clamps require use of MAGic Clamp Platform.



H1000-MR-CMB

MC-24™ Touch High Speed Microcentrifuge



- Full Color Touch Screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips (lid included)
- Compact design, 9 x 12 x 7.75 in.
- Unique air-flow system, keeps samples cool
- Optional rotors for 5ml and cryo tubes

OVER 30% OFF! SAVE OVER \$900...



C2417	MC-24™ Touch High Speed Microcentrifuge	\$2,699 \$1,795
C2417-R5	Optional 8 x 5ml rotor for MC-24™ Touch	\$1,150
C2417-R5-ADP2	Cryovial adapter for C2417-R5 rotor, 8/pk	\$144

BLACK FRIDAY DEALS



MULTI-THERM™
COOL - HEAT - SHAKE *Touch*

- Large touch screen for simple and easy programming and viewing of status
- Advanced programming features include adjustable ramp rates, preheat and clockwise/counterclockwise mixing
- Multiple block options for plates and tubes up to 50ml
- Broad temperature and speed range to meet a variety of applications

H51000-HCT Multi-Therm Touch Shaker with heating/cooling ~~\$4,421~~ **\$3,921**
H5100-CMB Combination block: 15 x 0.5ml & 20 x 1.5ml **\$394**
H5100-HL Heated lid for Multi-Therm Touch **\$844**
 Blocks sold separately. Additional blocks available. Ask for details.



**LIMITED
TIME
ONLY!
ORDER
TODAY!**

SMARTREADER™ 96

- Stand alone operation - no computer needed
- User- friendly, touchscreen programming
- Integrated plate shaker
- Supplied with 405, 450 , 492 and 630nm filters



MR9600 SmartReader 96 Microplate Reader ~~\$7,103~~ **\$6,103**
MR9600-TP Thermal Printer for MR9600 **\$373**
MR9600-PA Extra paper for printer, pk of 4 **\$36**



8 CHANNEL



12 CHANNEL

SMARTWASHER™ 96

- Precision needle position control
- Auto-rinse to eliminate cross contamination
- Liquid level sensors in bottle
- Customize and store up to 100 wash protocols



MW9600 Microplate Washer without manifolds ~~\$5,473~~ **\$4,473**
MW9600-12CM 12 channel washer manifold **\$389** **FREE!**
MW9600-8CM 8 channel washer manifolds **\$353** **FREE!**

Free offers require purchase of designated items at listed price. **Offers expire 12/05/2025.** Not valid on previous purchases. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code **TURKEY25**

BENCHMARK SCIENTIFIC PHONE (908)769-5555 EMAIL: INFO@BENCHMARKSCIENTIFIC.COM