

MC-24 Series High Speed Microcentrifuges

C2417



MC-24 Touch Ambient

- Full Color Touch Screen
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes and strips
- Unique air-flow system, keeps samples cool

Ambient
C2417 **\$2,693**

C2417-R



REFRIGERATED!

MC-24R Refrigerated

- Touch operation
- Powerful refrigeration, as low as -20°C
- COMBI-Rotor™ accepts tubes and strips
- 16,800 xg!
- Status indicator light

Refrigerated
C2417-R **ONLY \$4,999!***
*promo code CEN26



Both models included with **Combi-Rotor™**

Accepts:
24 x 1.5/2.0ml
2x PCR strips
(Lid included)

Hermle-Benchmark Universal Centrifuges

(Additional models available, see pg.2)

Z327-K



Z327 Series Bundles

- **Tissue culture:** Up to 8 x 50ml and 16 x 15ml
- **Blood tube:** Up to 28 x 8-10ml and 32 x 1.8-7ml



Tissue Culture Pkg.



Blood Tube Pkg.

5,000 rpm / 4,136 xg

Ambient
Z327-BND
(Tissue Culture)
Z327-BND-BD
(Blood tube) **\$9,700**

Refrigerated
Z327-K-BND
(Tissue Culture)
Z327-K-BND-BD
(Blood tube) **\$13,000**

Z446-K




Z446 Series Bundles

- **Tissue culture:** Up to 28 x 8-10ml and 32 x 1.8-7ml
- **Microplate:** Up to 20 std. plates and 4 x Deepwell



Tissue Culture Pkg.



Microplate Pkg.

Ambient:
4,000 rpm / 3,434 xg

Refrigerated:
4,500 rpm / 4,346 xg



Ambient
Z446-BND
(Tissue Culture)
Z446-MPBND
(Microplate) **\$12,700**

Refrigerated
Z446-K-BND
(Tissue Culture)
Z446-K-MPBND
(Microplate) **\$16,700**

Need help choosing the right centrifuge?
Try our centrifuge selector app!



Hermle-Benchmark Compact Centrifuges

Z206-A



Z206-A Compact

- 6 x 50ml rotor



6,000 rpm / 4,427 xg

Z206-BND: **\$2,950**

Z307



Z307 Deepwell Plate

- 2 x Deepwell Plate rotor



4,500 rpm / 2,716 xg

Z307-BND-MP: **\$6,400**

PlateFuge Microplate Microcentrifuge

C2000



- Quickly spins down droplets in PCR plates
- Unique rotor design prevents sample spillage (patent pending)
- 50% smaller than traditional centrifuges
- Accepts all popular PCR plates (skirted, semi, non) and most microtiter plates

Unique swing out rotor



Start/End



During run

C2000: **\$894**

Sprint 6H PLUS Clinical Centrifuge

C3303-6HP



- Includes 6-place swing-out rotor for tubes up to 15ml
- 4,000 rpm / 2,504 xg
- Perfect for clinical applications
- Adjustable accel and decel for gentle separation
- Optional 12-place fixed angle rotor



Includes 6-place Swing-out rotor!

C3303-6HP: **\$1,898**

Additional Hermle-Benchmark Bundles

Z496-K Floor Standing

w/ 40 x 50ml / 96 x 15ml SO rotor



SIEVA III Filtration Centrifuge

w/ perforated or non perforated basket



Z366 Series Universal Centrifuges

w/ 16 x 50ml / 40 x 15ml SO rotor



Z36-HK Super Speed Centrifuge

w/ 3 rotor options up to 65,395 xg!

