

Offers valid until Dec. 31st 2026

Hermle-Benchmark Universal Centrifuges

(Additional models available, see pg.2)

Z327-K



Z327 Series Bundles

- **Tissue culture:** Up to 8 x 50ml and 16 x 15ml
- **Blood tube:** Up to 28 x 8-10ml and 32 x 1.8-7ml



Tissue Culture Pkg.



Blood Tube Pkg.

5,000 rpm / 4,136 xg

Ambient
Z327-BND
(Tissue Culture)
Z327-BND-BD
(Blood tube)

\$9,700

Refrigerated
Z327-K-BND
(Tissue Culture)
Z327-K-BND-BD
(Blood tube)

\$13,500
BenchPACK™

Z366-K



Z366 Series Bundles

- **Tissue culture:** Up to 16 x 50ml and 40 x 15ml



Ambient:
4,500 rpm / 3,735 xg

Ambient
Z366-BND

\$10,300

Refrigerated
Z366-K-BND

\$14,700
BenchPACK™

Z446-K



Z446 Series Bundles

- **Tissue culture:** Up to 28 x 50ml and 68 x 15ml
- **Microplate:** Up to 20 std. plates and 4 x Deepwell



Tissue Culture Pkg.



Microplate Pkg.

Ambient:
4,000 rpm / 3,434 xg

Refrigerated:
4,500 rpm / 4,346 xg

Ambient
Z446-BND
(Tissue Culture)
Z446-MPBND
(Microplate)

\$12,700

Refrigerated
Z446-K-BND
(Tissue Culture)
Z446-K-MPBND
(Microplate)

\$16,990
BenchPACK™



FREE

BenchPACK™

with refrigerated universal and super speed centrifuge bundles!
(Z327-K, Z366-K, Z446-K, Z496-K, Z36-HK bundles)



INCLUDES:

myFuge™ 12 Mini Centrifuge
BenchMixer™ V2 Vortex Mixer
Roto-Mini™ Rotating Mixer
MiniMag™ Mini Magnetic Stirrer
Accuris™ Mini Balance

Value \$1,750!

Need help choosing the right centrifuge?
Try our centrifuge selector app!



Offers valid until Dec. 31st 2026

PlateFuge 2D Digital Plate Centrifuge

C2000



- Digital control, adjustable rpm and time
- Swing out rotor prevents sample spillage (no plate seal required)
- Personal sized, 9 x 12" footprint
- Accepts all popular PCR plates (skirted, semi, non) and most microtiter plates

2,800 rpm / 728 xg

Unique swing out rotor



Start/End



During run

C2006: **\$895**

MC-24 Series High Speed Microcentrifuges

C2417

C2417-R



REFRIGERATED

- Full Color Touch Screen
- 16,800 xg
- COMBI-Rotor™ accepts tubes and strips
- Powerful refrigeration, as low as -20°C (C2417-R)
- Status indicator light (C2417-R)
- Additional rotor options



Ambient
C2417

ONLY \$1,990!*

Refrigerated
C2417-R

ONLY \$4,999!*

*promo code CEN26



Both models included with **Combi-Rotor™**

Accepts:
24 x 1.5/2.0ml
2x PCR strips
(Lid included)

Hermle-Benchmark Compact Centrifuges

Z206-A



Z287-A Micro Bundle

- 2 x microplate and 18 x 1.5/2.0ml



3,800 rpm / 3,717 xg



14,000 rpm / 15,994 xg

Z287-BND-MB: **\$5,500**

Z307



Z307 Deepwell Plate

- 2 x Deepwell Plate rotor



4,500 rpm / 2,716 xg

Z307-BND-MP: **\$6,400**

Additional Hermle-Benchmark Bundles

Z496-K Floor Standing

w/ 40 x 50ml / 96 x 15ml SO rotor



BenchPACK™

SIEVA III Filtration Centrifuge

w/ perforated or non perforated basket



Z216 Series Microcentrifuges

w/ 2 rotor options up to 21,631 xg!



Z36-HK Super Speed Centrifuge

w/ 3 rotor options up to 65,395 xg!



BenchPACK™