

Offers valid until Dec. 31st 2026

## Hermle-Benchmark Universal Centrifuges

(Additional models available, see pg.2)

Z327-K



### Z327 Series Bundles

- **Tissue culture:** Up to 8 x 50ml and 16 x 15ml
- **Blood tube:** Up to 28 x 8-10ml and 32 x 1.8-7ml



Tissue Culture Pkg.



Blood Tube Pkg.

5,000 rpm / 4,136 xg

**Ambient**  
Z327-BND  
(Tissue Culture)  
Z327-BND-BD  
(Blood tube)

**\$9,700**

**Refrigerated**  
Z327-K-BND  
(Tissue Culture)  
Z327-K-BND-BD  
(Blood tube)

**\$13,500**

**BenchPACK™**

Z366-K



### Z366 Series Bundles

- **Tissue culture:** Up to 16 x 50ml and 40 x 15ml



Ambient:  
4,500 rpm / 3,735 xg

**Ambient**  
Z366-BND

**\$10,300**

**Refrigerated**  
Z366-K-BND

**\$14,700**

**BenchPACK™**

Z446-K




### Z446 Series Bundles

- **Tissue culture:** Up to 28 x 8-10ml and 32 x 1.8-7ml
- **Microplate:** Up to 20 std. plates and 4 x Deepwell



Tissue Culture Pkg.



Microplate Pkg.

Ambient:  
4,000 rpm / 3,434 xg

Refrigerated:  
4,500 rpm / 4,346 xg

**Ambient**  
Z446-BND  
(Tissue Culture)  
Z446-MPBND  
(Microplate)

**\$12,700**

**Refrigerated**  
Z446-K-BND  
(Tissue Culture)  
Z446-K-MPBND  
(Microplate)

**\$16,990**

**BenchPACK™**



**FREE**  
**BenchPACK™**

with refrigerated universal and super speed centrifuge bundles!  
(Z327-K, Z366-K, Z446-K, Z496-K, Z36-HK bundles)



### INCLUDES:

- myFuge™ 12 Mini Centrifuge
- BenchMixer™ V2 Vortex Mixer
- Roto-Mini™ Rotating Mixer
- MiniMag™ Mini Magnetic Stirrer
- Accuris™ Mini Balance

Value **\$1,750!**

Need help choosing the right centrifuge?  
Try our centrifuge selector app!



Offers valid until Dec. 31st 2026

## PlateFuge 2D Digital Plate Centrifuge

C2000



- Digital control, adjustable rpm and time
- Swing out rotor prevents sample spillage (no plate seal required)
- Personal sized, 9 x 12" footprint
- Accepts all popular PCR plates (skirted, semi, non) and most microtiter plates

2,800 rpm / 728 xg

Unique swing out rotor



Start/End



During run

C2006: **\$895**

## MC-24 Series High Speed Microcentrifuges

C2417

C2417-R



**REFRIGERATED**

- Full Color Touch Screen
- 16,800 xg
- COMBI-Rotor™ accepts tubes and strips
- Powerful refrigeration, as low as -20°C (C2417-R)
- Status indicator light (C2417-R)
- Additional rotor options



Ambient  
C2417

**ONLY \$1,990!\***

Refrigerated  
C2417-R

**ONLY \$4,999!\***

\*promo code CEN26



Both models included with **Combi-Rotor™**

Accepts:  
24 x 1.5/2.0ml  
2x PCR strips  
(Lid included)

## Hermle-Benchmark Compact Centrifuges

Z206-A

### Z287-A Micro Bundle

- 2 x microplate and 18 x 1.5/2.0ml



3,800 rpm / 3,717 xg




14,000 rpm / 15,994 xg

Z287-BND-MB: **\$5,500**

Z307

### Z307 Deepwell Plate

- 2 x Deepwell Plate rotor



4,500 rpm / 2,716 xg

Z307-BND-MP: **\$6,400**

## Additional Hermle-Benchmark Bundles

### Z496-K Floor Standing

w/ 40 x 50ml / 96 x 15ml SO rotor



**BenchPACK™**

### SIEVA III Filtration Centrifuge

w/ perforated or non perforated basket



### Z216 Series Microcentrifuges

w/ 2 rotor options up to 21,631 xg!



### Z36-HK Super Speed Centrifuge

w/ 3 rotor options up to 65,395 xg!



**BenchPACK™**