Benchmark Scientific

Fall into SAVINGS

BLACK FRIDAY SALE

Offers expire 12/04/2026
Use promo code: TURKEY26



PlateFuge™ 2D
Digital MicroPlate Centrifuge

FREE MYFUGE MINI With Purchase!
$347 VALUE

- Digital control, adjustable rpm and time
- Swing out rotor prevents sample spillage (no plate seal required)
- Personal sized, 9 x 12" footprint
- Accepts all popular types of PCR plates (skirted, semi, and non) and most microtiter plates

C2006	*PlateFuge™ 2D Digital MicroPlate MicroCentrifuge*	**$939.00**
C1008	*myFuge™ Mini Centrifuge*	**FREE**



C1008



ACCURIS INSTRUMENTS COMPACT BALANCE

- User friendly interface
- Easy in-lab calibration feature
- RS232 port for data output
- 500g capacity / 0.01g readability

W3300-500	*Accuris™ Compact Balance*	**$412.92**
C1008	*myFuge™ Mini Centrifuge*	**FREE**

BUY ONE GET ONE *FREE!*

Roto-Mini™ Plus

- Included tube holders:
  - 11 x 1.5/0.2ml (pack of 2)
  - 5 x 10-15ml (pack of 2)
  - 5 x 5-7ml (pack of 2)
  - 6 x 50ml (pack of 1)
- Variable speed control, 5 to 70 rpm
- Advanced shaking/mixing modes: Mix (oscillating), Pause, Rock



R2024	*Roto-Mini™ Plus Rotator, Variable Speed*	**$541.32**

FOR TWO!



Reversible Head

BenchMixer™ V2

- Reversible top platform, cup or flat-dimpled surface
- Powerful motor for INSTANT vortexing
- Continuous or "touch" operation
- Q-Drive™ dynamic balancing system
- Variable speed from 200 to 3200 rpm
- Wide variety of accessories, for tubes up to 50ml

BUY ONE GET ONE *FREE!*

BV1003	*BenchMixer V2, Vortex Mixer with flip top cup head*	**$342.08**

BENCHMARK SCIENTIFIC PHONE (908)769-5555 | EMAIL: INFO@BENCHMARKSCIENTIFIC.COM

SALE

BLACK FRIDAY SALE



myFUGE™12 Plus


- COMBI-Rotor for tubes **and** strips
- Digital control up to 7,000 rpm / 2,739 xg!
- Twice the capacity of traditional mini centrifuges
- Near silent operation
- Timer up to 99 mins.

C2212 *MyFuge™12 Plus Digital Mini Centrifuge* **$708.72**
FOR TWO!

BUY ONE GET ONE FREE!

FREE 500g of Agarose!





Includes two casting sets

MY GEL™ mini

- Complete system, includes gel box, power supply and casting sets
- Quick and easy gel casting, no tape, clamps or spacers required
- Compact design, footprint: 7 x 5.1 in.

E1101 *myGel™ Mini Electrophoresis System* **$691.86**
A1705 *Benchmark Agarose LE, 500g* **FREE**



A1705



BEADBUG™3 MICROTUBE HOMOGENIZER



- Extremely powerful mixing for lysis, grinding or homogenization
- Faster and more effective than tissue grinders
- Extremely compact, less than 7" wide
- Eliminates the cross contamination associated with handheld homogenizers
- Full range of beads in prefilled tubes and bulk available

D1030 *BeadBug™3 Microtube homogenizer, 115V* **$1,022.95**
FOR TWO!

Buy One Get One Free!



Standard Glass Beads (Acid Washed & Heat Treated) | Triple-Pure™ (Molecular Biology Grade) Zirconium Beads | Stainless Steel Beads (Alcohol Washed & Heat Treated)

A variety of bead sizes to support different sample types and applications are available.
For customized recommendations or specific application support, contact our team directly at **info@benchmarkscientific.com**.

BENCHMARK SCIENTIFIC PHONE (908)769-5555 | EMAIL: INFO@BENCHMARKSCIENTIFIC.COM

Benchmark Scientific

Free offers require purchase of designated items at listed price. ***Offers expire 12/04/2026.*** Not valid on previous purchases. Special price applies only to item specified. Use promo code: ***TURKEY26***

BLACK FRIDAY

30% SALE



SMARTWASHER™ 96

- Precision control of aspiration needle position
  - Accurate to 0.1mm
  - Minimal residual liquid, ≤1µl/well
- Auto rinse function eliminates cross contamination
- Up to 100 wash protocols can be customized and stored
- Liquid level sensors in wash and waste bottles
- Wash protocols can be transferred via flash drive

More Great Deals For Your Lab This Season!

MW9600	*SmartWasher™ 96 Plate Washer*	~~$5609.72~~	**$3926.80**
MW9600-12CM	*12 Channel Manifold for MW9600*	**$404.48**	
MW9600-8CM	*8 Channel Manifold for MW9600*	**$367.08**	

Over 30% OFF!

INCU-SHAKER™ Mini HC

ECO FRIENDLY



- Heating **and Cooling**, -8°C from ambient to 50°C
- Extremely compact, just 11" wide
- TOUCH Screen control over time, speed and temperature
- Eco-friendly **HFC free cooling**
- Instantly exchange flask clamps with MAGic Clamp™ accessories (PATENTED)

H1001-M-HC	*Incu-Shaker™ Mini HC heating and cooling shaking incubator*	~~$3,996.00~~	**$ 2797.20**
H1000-MR	*MAGic Clamp universal platform, (9.5 x 11.5")*	**$326.70**	
H1000-MR-CMB	*MAGic Clamp™ COMBO, for 125ml, 250ml, 500ml flasks*	**$58.64**	

MC-24R™ Refrigerated Microcentrifuge





COMBI-Rotor (Included)

- Included COMBI-Rotor™ accepts tubes **and** strips (lid included)
- Powerful refrigeration, as low as -20°C
- Touch Operation with speed control to 16,800 xg
- Status indicator light

C2417-R	*MC-24R™ Refrigerated High Speed Microcentrifuge*	~~$6,087.87~~	**$4261.51**

ACCURIS™ LIFE SCIENCE REAGENTS — qMAX™ Green qPCR Mix With Blue Tracking Dye

BOGO

- Superior sensitivity and fast cycling with exceptional results
- Includes Accuris Hot Start Taq for greater specificity and accuracy
- Ideal for low copy number templates
- 2X master mix with inert blue tracking dye



PR2000-N-100	*Accuris qMAX Green qPCR, No ROX, 100 rxns*	**$83.67**	**FOR TWO!**
PR2000-L-100	*Accuris qMAX Green qPCR, Low ROX, 100 rxns*	**$83.67**	**FOR TWO!**
PR2000-H-100	*Accuris qMAX Green qPCR, High ROX, 100 rxns*	**$83.67**	**FOR TWO!**

BENCHMARK SCIENTIFIC PHONE (908)769-5555 | EMAIL: INFO@BENCHMARKSCIENTIFIC.COM